

**The value of research
protected by strict
production management**

BlackCheckTM

***MG(Mycoplasma gallisepticum)* Differential Identification Kit**

A multi-gene diagnostic kit that rapidly and accurately identifies four specific genes of *Mycoplasma gallisepticum* (3-vaccine strains: ts-11, 6/85, F & 1-field strain: MG), the causative agent of respiratory diseases in poultry

www.blackpcr.com

MG (*Mycoplasma gallisepticum*) Differential Identification Kit

Mycoplasma gallisepticum is commonly involved in the polymicrobial "chronic respiratory disease" in broiler chickens, leading to increased condemnations in the processing plant.

In layers and breeders, it is usually subclinical, but causes a reduction in the number of eggs laid per hen over the production cycle. In chickens, this disease not only causes respiratory symptoms but also leads to reduced feed efficiency, lower weight gain, decreased egg production, and hatchability.

MG (*Mycoplasma gallisepticum*) Differential Identification Kit is a Detection type product that simultaneously rapidly and accurately identifies four causative genes of *Mycoplasma gallicepticum* that cause respiratory diseases in poultry.

In addition to primers that specifically bind to *Mycoplasma gallisepticum* causative genes (a total of 4 types of vaccine strain and field strain), it contains all the components necessary for PCR execution, so it is possible to confirm the presence or absence of infection and identification results simply by adding template DNA and D.W.

Product Data






Figure 1. MG (*Mycoplasma gallicepticum*) Differential Identification Kit repeatability (three repetitions) test results

Order Details

Cat No.	Product Name	Specification	
BMGK-48T	MG (<i>Mycoplasma gallisepticum</i>) Differential Identification Kit	48 test (8 strip tube x 6 ea) 15 µl/tube	20µl reaction

Distributor :  [41516] 33, Yutongdanji-ro 24-gil, Buk-gu, Daegu, Republic of Korea
 Tel: 053 383 2989 Fax: 053 383 2986 E-mail: kc1111@naver.com

Manufacturer :  [08592] 917, 9F, 14, Gasan digital 2-ro, Geumcheon-gu, Seoul, Republic of Korea
 Tel: 02 2061 8033 Fax: 02 2061 8034 E-mail: sales2@forbio.co.kr
 web: www.blackpcr.com

