

The value of research
protected by strict
production management

BlackCheckTM

NDV Duplex RT-qPCR Kit

**A rapid and accurate Newcastle Disease Virus(NDV) diagnostic kit
capable of simultaneously detecting pathogenic
and non-pathogenic strains**

www.blackpcr.com

NDV Duplex RT-qPCR Kit


Newcastle disease (ND) in poultry is caused by the Newcastle disease virus (NDV).

The disease is highly contagious, can cause respiratory, gastrointestinal and neurological symptoms when infected, and has a mortality rate of nearly 100%, causing significant economic damage to the poultry industry.

NDV is classified into velogenic, mesogenic, lentogenic, and apathogenic depending on the WOAHP considers infection by NDV of mesotype or higher as an outbreak of Newcastle disease.

The NDV Duplex RT-qPCR Kit simultaneously detects pathogenic and non-pathogenic properties, and can quickly and accurately diagnose the occurrence of Newcastle disease.

Product Data



Order Details

Cat No.	Product Name	Specification	
BNDVK-48T	NDV Duplex RT-qPCR Kit	48 test (8 strip tube x 6 ea) 15 µl/tube	20µl reaction

Distributor :



[41516] 33, Yutongdanji-ro 24-gil, Buk-gu, Daegu, Republic of Korea
Tel: 053 383 2989 Fax: 053 383 2986 E-mail: ventech2989@naver.com

Manufacturer :



08592] 917, 9F, 14, Gasan digital 2-ro, Geumcheon-gu,
Seoul, Republic of Korea
Tel: 02 2061 8033 Fax: 02 2061 8034 E-mail: sales2@forbio.co.kr

ISO 13485:2016 Design and Development, Manufacturing, Distribution and Sales, and Service of In vitro Diagnostic Reagents (Molecular Diagnostic Devices)

